

DECUS

PROGRAM LIBRARY

DECUS NO.	12-69
TITLE	AN ON-LINE FOCAL-12 PROGRAM FOR AUTO-ANALYZERS
AUTHOR	Mack W. Overton, Jr., Larry L. Alber and Dr. Donald E. Smith
COMPANY	U. S. Food and Drug Administration Chicago, Illinois
DATE	March 10, 1972
SOURCE LANGUAGE	FOCAL-12

AN ON-LINE FOCAL-12 PROGRAM FOR AUTO-ANALYZERS

DECUS Program Library Write-up

DECUS NO. 12-69

A program and hardware interface have been developed to acquire and reduce data from an auto-analyzer. A system diagram is presented in Fig. 1. A PDP-12A (8K core) was used for the project. The interfacing circuitry is shown in Fig. 2 and is based on the use of a retransmitting slidewire mounted on the pen cable pulley of the spectrophotometer recorder. The signal taken off the center tap of the slidewire is then conditioned as shown and fed to the A/D converter inputs of the PDP-12A.

The program is written in a modified version of FOCAL-12 (1) and utilizes FOCAL-12 library capabilities for program and user storage. The program is contained in four segments on magnetic tape unit 0. The first segment is loaded via teletype command and the remaining segments are loaded under program control. Two files on magnetic tape unit 1 provide storage for user input and a maximum of 70 readings (samples and standards) taken from the auto-analyzer. Although storage for a maximum of 70 readings is provided, the first two and the final reading are ignored by the calculation routines. Therefore, in order to utilize fully the storage provided, some program modification would have to be made.

%ONITA1, the first segment, is an initialization routine which is loaded via the library go command issued from the teletype by the user. It is during initialization that basic information concerning the run is supplied by the user. This is accomplished through a user interactive scope-teletype routine.

%ONITA2, the second segment, is loaded under program control upon completion of %ONITA1 and consists of the realtime sample, display and peak hunt routines. A sense switch (sns 1) on the computer console controls the exit and entry into the peak hunt routine. Only the X and Y coordinates for peaks, initial baseline and final baseline are stored. Raw data are not stored.

When sampling is complete, %ONITA3 automatically calculates and outputs the results. At this point the program determines if it found the specified number of peaks, and the user is asked to verify peak identity. If an error is detected %ONITA4 is loaded and the user must code the absorbance values as the computer outputs them by typing S (standard), U(unknown), C(composite), or X (deletion). It should be noted that while %ONITA3 uses an average standard value, %ONITA4 calculates using individual standards (i.e., the most recent standard). Baseline correction is accomplished by fitting a point to a straight line constructed between the initial and final baseline points. This value is subtracted from the corresponding absorbance prior to calculation.

¹ FOCAL-12 with overlay to permit reading of sense switches on the PDP-12A console. Courtesy of Mr. Dale L. Lewellyn, Software Consultant, Digital Equipment Corporation, 1350 Frontage Road, Northbrook, Illinois.

A versatile inexpensive system utilizing the PDP-12A has been developed for on-line data acquisition and reduction for auto-analyzers. While admittedly deficient in two areas, no provision for raw data storage or time sharing, it should find use in projects of limited duration which prohibit development of a more sophisticated but expensive system. It also serves admirably as an interim measure during the development of a more sophisticated system.

NOTE

The use of equipment described in this paper does not constitute endorsement by the Food and Drug Administration.




FIG. 1. SYSTEM SCHEMATIC




FIG. 2. Interface circuit constructed with Burr-Brown Model 603 operational amplifier manifold and plug in components, Burr-Brown Research Corp., Tucson, Arizona 85706

*L L,%ONITA1,Ø

*W A

C FOCAL-12

Ø1.Ø2 L M, 1, ITA, 1;L O, F1, F, ITA, 1:F I=Ø,84;S F1(I)=Ø
Ø1.Ø4 L M, 1, ITA2, 1;L O, F3, F, ITA2, 1;F I=Ø,84;S F3(I)=Ø
Ø1.Ø6 O C; O S
Ø1.Ø7 A !, "CH NO", F1(Ø)
Ø1.Ø9 A !, "NO OF ITAS, EXCLUDE COMP", F3(7Ø)
Ø1.1Ø A !, "NO PEAKS, SPLS+STDS+COMP", F3(71)
Ø1.11 A !, "COMP CUP NO. ", F3(72)
Ø1.12 O C;DO 2;A !, "STD CONC", F3(8Ø)
Ø1.13 A !, "DIL FAC", F3(81)
Ø1.14 A !, "DECLARED IN MG OR GR?", MG
Ø1.15 I (MG-ØMG)1.17, 1.16, 1.14
Ø1.16 S F3(82)=1;GOTO 1.18
Ø1.17 S F3(82)=1/64.8
Ø1.18 A !, "AMT DECLARED", F3(74)
Ø1.19 A !, "WT OF COMP", F3(83)
Ø1.2Ø A !, "AVG TAB WT", F3(84)
Ø1.21 S F3(73)=F3(8Ø)*F3(81)*F3(82)
Ø1.22 S F3(75)=F3(83)/F3(84)
Ø1.23 L C, F1;L C, F2;L C, F3
Ø1.24 L G, %ONITA2, Ø
Ø1.3Ø T "ERROR "
Ø1.31 Q

Ø2.Ø3 A !, "NO. OF STDS", F1(7Ø)
Ø2.Ø4 T !, "STD CUP NO"
Ø2.Ø5 F I=1, F1(7Ø);DO 3.Ø2
Ø2.Ø6 A !, "STD DATA OK, Y OR N?", Y
Ø2.Ø7 I (Y-ØY)2.Ø8, 2.Ø9, 2.Ø8
Ø2.Ø8 O C;GOTO 2.Ø3
Ø2.Ø9 O C

Ø3.Ø2 T !, I, " ";A F1(I+7Ø)
*

*L L,%ONITA2,Ø

*W A

C FOCAL-12

```
Ø1.Ø1 A "T",T;O C
Ø1.Ø2 L O,F1,F,ITA,1;L O,F3,F,ITA2,1
Ø1.Ø3 O C;O I,T;S J=Ø;S K=Ø;S P=Ø;S M=Ø;S F3(81)=-1
Ø1.Ø4 DO 2;DO 2
Ø1.Ø5 I (FX(1))1.4Ø,1.Ø4,1.Ø6
Ø1.Ø6 S F3(77)=<F3(83)+F3(84)>/2;DO 1.Ø3
Ø1.Ø7 I (FX(1)-1)1.3Ø;DO 2
Ø1.Ø8 I (F3(84)-F3(83))1.1Ø,1.1Ø
Ø1.Ø9 S P=P+1;I (P-5)1.11,1.12,1.4Ø
Ø1.1Ø S P=Ø
Ø1.11 GOTO 1.Ø7
Ø1.12 DO 1.22;S F3(81)=F3(84)
Ø1.14 I (FX(1)-1)1.3Ø;DO 2
Ø1.16 I (F3(84)-F3(83))1.17,1.18,1.18
Ø1.17 S M=M+1;I (M-5)1.19,1.2Ø,1.4Ø
Ø1.18 S M=Ø
Ø1.19 GOTO 1.14
Ø1.2Ø S F3(76)=F3(76)+1
Ø1.21 S J=J+1;S F1(J)=F3(82);S F3(J)=F3(81)
Ø1.22 S P=Ø;S M=Ø
Ø1.23 GOTO 1.Ø7
Ø1.3Ø S F3(78)=K;S F3(79)=<F3(83)+F3(84)>/2
Ø1.33 L C,F1;L C,F3;L G,%ONITA3,Ø
Ø1.4Ø T "ERR";Q

Ø2.Ø1 S K=K+(1Ø*T)
Ø2.Ø2 S F3(83)=F3(84)
Ø2.Ø3 S Q=Ø;F I=1,1Ø;S Q=Q+FADC(F1(Ø));O I
Ø2.Ø4 S F3(84)=Q/1Ø
Ø2.Ø5 I (F3(84)-F3(81))2.Ø6;S F3(81)=F3(84);S F3(82)=K
Ø2.Ø6 S O=FDIS(K/(1ØØØØ*T),F3(84))
*
```


*L L,%ONITA3,Ø

*W A

C FOCAL-12

Ø1.Ø3 O T
Ø1.Ø4 L O, F1, F, ITA, 1
Ø1.Ø5 L O, F3, F, ITA2, 1
Ø1.Ø9 S $M = \langle F3(79) - F3(77) \rangle / F3(78); S ST = \emptyset$
Ø1.1Ø F $I = 73, 7Ø + F1(7Ø) - 1; S ST = ST + F3(F1(I)) - (\langle M * F1(F1(I)) \rangle + F3(77))$
Ø1.12 S $ST = ST / \langle F1(7Ø) - 3 \rangle$
Ø1.14 T !, "SECONDS NET PH CODE FOUND % DECLARED", !, !
Ø1.16 S $I = \emptyset; S SN = 71; S QA = \emptyset$
Ø1.17 S $I = I + 1; I (F3(76) - I) 1.22$
Ø1.18 S $SA = F3(I) - (\langle M * F1(I) \rangle + F3(77))$
Ø1.19 T !, %4.Ø, F1(I), %8.Ø3, SA
Ø1.2Ø I $(F1(SN) - I) 1.4Ø, 2.Ø1$
Ø1.21 I $(I - F3(72)) 2.1Ø, 2.2Ø, 2.1Ø$
Ø1.22 T !, "AVG(", %2.Ø, F3(7Ø), " " , %8.Ø3, QA/F3(7Ø)
Ø1.23 T $(QA / F3(7Ø)) * (1ØØ / F3(74))$
Ø1.24 A !, "RECALC, Y OR N?", Y
Ø1.25 I $(Y - ØY) 1.26, 1.41, 1.24$
Ø1.26 T !, "DECL-", F3(74)
Ø1.27 T !, "CONC S X DIL X UNITS FAC =", F3(73)
Ø1.28 T !, "COMP/AVG TAB =", F3(75)
Ø1.29 T !, "BL1 =", F3(77), !, "BL2 =", F3(79), !, "TIME =", F3(78)
Ø1.3Ø A !, "SPL NO", Q
Ø1.31 A !, "DATE", Q
Ø1.32 A !, "INITIALS", Q
Ø1.33 A !, "PRODUCT", Q
Ø1.34 Q
Ø1.4Ø T !, "ERROR"; Q
Ø1.41 L G, %ONITA4, O

Ø2.Ø1 T " S"; S SN=SN+1; GOTO 1.17

Ø2.1Ø T " U"

Ø2.12 S $Q = \langle SA / ST \rangle * F3(73)$

Ø2.13 S $QA = QA + Q$

Ø2.14 GOTO 2.23

Ø2.2Ø T " C"

Ø2.21 DO 2.12

Ø2.22 S $Q = Q / F3(75)$

Ø2.23 T $Q, \langle Q / F3(74) \rangle * 1ØØ$

Ø2.24 GOTO 1.17

*

PEAK DETECTION

Y-COORDINATE

PEAK SENSED

AVG = BL1

SNS I HI

BASELINE CORRECTION

AVG = BL2

SNS I LO

END OF PEAK

X-COORDINATE

NOTE: Absorbance is converted to volts DC by the retransmitting slidewire connected to recorder pen slidewire pulley, and DC volts are in fact used for calculations.

DEK 50120 10 INCR

A
B
S
O
R
B
A
N
C
E
*

T I M E

DATA STORAGE SCHEMATIC

FILE F1		FILE F2	
LOCATION	CONTENT	LOCATION	CONTENT
∅	a/d ch. no.	∅	(vacant)
1	x-coord. peak1	1	y-coord. peak 1
2	x-coord. peak2	2	y-coord. peak2
3	x-coord. peak3	3	y-coord. peak3
4	x-coord. peak4	4	y-coord. peak4
5	x-coord. peak5	5	y-coord. peak5
.	etc.	.	etc.
.	etc.	.	etc.
.		.	
7∅	# of stds.	7∅	# ITAS (tablets)
71	first std cup #	71	# cups (total)
72	2nd std cup #	72	comp. cup #
73	3rd std cup #	73	(std conc) (dil fac) (1/units)
74	etc.	74	amt. declared
75	etc.	75	wt. of comp./avg. tab. wt.
76	etc.	76	# of peaks found by program
77	etc.	77	y-coord. BL1 (x=∅ by defn.)
78	"	78	x-coord. BL2.
79	"	79	y-coord. BL2.
8∅	"	8∅	Std conc
81	"	81*	std conc and max. y value
82	"	82*	units and x-coord. of y max.
83	"	83*	comp. wt. and current y
84	"	84*	avg. tab. wt. and current x

*(See notes below.)

NOTES

- Files are floating point format files. Each file location actually occupies three core locations.
- File locations 81 - 84 are multi-purpose utility storage locations. The first usage indicates contents and purpose for %ONITA1 and the second usage indicates contents and purpose for %ONITA2.

SECONDS	NET PH	CODE	FOUND	% DECLARED
339	∅.565	S		
459	∅.581	S		
573	∅.589	S		
693	∅.562	U	47.966	95.931
819	∅.513	U	43.806	87.612
936	∅.567	U	48.369	96.738
1053	∅.524	U	44.709	89.418
1176	∅.565	U	45.205	96.410
1290	∅.574	S		
1413	∅.551	U	47.040	94.080
1539	∅.571	U	48.751	97.503
1653	∅.535	U	45.642	91.283
1773	∅.549	U	46.886	93.771
1893	∅.567	U	48.396	96.793
2019	∅.589	S		
2136	∅.561	U	47.916	95.831
2253	∅.554	U	47.308	94.616
2376	∅.523	U	44.649	89.298
2493	∅.560	U	47.778	95.555
2613	∅.548	U	46.753	93.506
2733	∅.590	S		
2853	∅.554	U	47.273	94.546
2976	∅.556	U	47.466	94.933
3090	∅.568	U	48.476	96.953
3216	∅.591	U	50.455	100.909
3330	∅.536	U	45.744	91.487
3453	∅.578	S		
3570	∅.564	U	48.132	96.263
3690	∅.621	U	53.045	106.090
3816	∅.565	U	48.252	96.503
3930	∅.518	U	44.258	88.516
4053	∅.550	U	46.953	93.906
4176	∅.585	S		
4287	∅.481	U	41.068	82.135
4410	∅.556	U	47.432	94.864
4533	∅.558	U	47.626	95.251
4653	∅.585	U	49.987	99.974
4770	∅.574	U	49.012	98.025
4890	∅.603	C	51.457	102.914
5013	∅.594	S		
5130	∅.597	S		
AVG(30)			47.312	94.623

RECALC, Y OR N?:N
DECL= 50.000
CONC S X DIL X UNITS FAC = 50.000
COMP/AVG TAB= 1.000
BL1=- 0.382
BL2=- 0.387
TIME= 5985.000
SPL NO:EXHIBIT
DATE:9MAR72
INITIALS:AUTHORS
PRODUCT:EXHIBIT *

DATA ACQUISITION AND CALCULATION

1. Data is acquired in groups of ten a/d samplings, each taken at the interval specified in the user's response to T during initialization. These ten readings are summed and the average is then considered as the y coordinate of a data point for defining curve shape. The x coordinate is time elapsed since setting SNS 1.
2. A peak is defined as five consecutive positive first derivatives followed after an unspecified interval of time by five consecutive negative first derivatives. The first derivative is obtained by subtracting the y coordinate of the previous data point from the y coordinate of the current data point. It should be noted that a zero first derivative will break the chain of five consecutive first derivatives.
3. A peak maximum is defined as the largest value obtained between the five consecutive positive and the five consecutive negative first derivatives.
4. Only the x and y coordinates for the peaks, initial baseline and final baseline are stored. All other data are lost.
5. The initial and final baseline y coordinates are defined as the average of the y coordinates of the last two data points taken prior to changing the position of SNS 1. (Note: x coordinate for BL1 = \emptyset).
6. Results are calculated from the following equations:

$$(a) \quad \text{Found} = \frac{\text{Abs}_u}{\text{Abs}_s} \times \text{Conc Std} \times \text{Dil} \times \frac{1}{\text{Units}}$$

Abs_u = Sample abs. - Baseline Correction

Abs_s = Standard abs. - Baseline Correction

Conc Std is expressed in mg/ml

Dil = Volume to which 1 tablet is diluted (mls.)

$1/\text{Units}$ = 1/1 for declared in mg.
1/64.8 for declared in grains (grs)

$$(b) \quad \% \text{ Declared} = \frac{\text{Found}}{\text{Declared}} \times 100\%$$

7. Baseline correction is accomplished by calculating the slope for a line between BL 1 and BL2. (See Eq. (c)).

$$(c) \quad \text{Slope (m)} = \frac{y_2 - y_1}{x_2 - x_1}$$

The slope is then applied in equation (d) to calculate the baseline correction for each peak.

$$(d) \quad y = mx + b; \text{ Where } y = \text{Baseline correction}$$
$$x = x \text{ coordinate of peak max.}$$
$$b = \text{BLI}$$

8. The value used for standard absorbance (Abs.) is obtained differently in %ONITA3 and %ONITA4. %ONITA3 determines an average standard value composed of all the standards excluding the first two standards and the final standard. %ONITA4 calculates using an individual standard value which must precede the unknowns or composite. See also next section.

SEQUENCING OF SAMPLES, STANDARDS AND COMPOSITE

1. Storage is provided for a maximum of 70 peaks, however this does not mean 70 unknowns. Careful study of the program and enclosed storage schematic will reveal two important conventions which must be adhered to in sequencing samples, standards and the composite in the auto-analyzer.

a. %ONITA3 requires a minimum of 4 standards, since the first two and the last standard value are not used for calculating.

b. %ONITA4 assumes at least one standard precedes any samples or the composite.

In view of the above, it becomes apparent that for the programs to work properly the run must contain at least 4 standards, one of which must be located in the first cup.

2. Provision is made for only one composite.

3. The above are minimum requirements and must be met for both programs to operate properly.

4. Sequencing of samples, standards and the composite beyond these requirements is left to the user's discretion.

Restriction: DIAL-MS modified so that if only one DF32 is present, Set Sense Switch \emptyset . For all other configurations, operate DIAL-MS as usual.

OPERATING INSTRUCTIONS

1. Mount DIAL-MS tape containing program on tape unit \emptyset and formatted LINCtape on unit 1.

2. Start DIAL-MS in the usual manner:

- A. LSW $\emptyset 7 \emptyset 1$ RSW 731 \emptyset
- B. Depress IO PRESET
- C. Depress DO
- D. Depress START $2\emptyset$ when the tape movements stop

3. When cursor appears on scope, type " \rightarrow LO FOCALSNS, $\emptyset \downarrow$ " to load modified version of FOCAL-12.

4. When TTY responds with an "*" type "L G,%ONITAI,Ø)" to load the first segment of the program.

5. Once %ONITAI is loaded the following display will begin. Answer all questions on TTY. Use decimal numbers. Avoid FOCAL-12 terminators during response. The program is formatted assuming the space bar will be used to terminate responses. Typing an "←" deletes everything to the left of the "←" within the response, up to the colon.

INITIAL DISPLAY

CH NO:	a/d channel
NO OF ITAS, EXCLUDE COMP:	# of tablets to be assayed
NO PEAKS, SPLS+STDS+COMP:	Total # of cups in run
COMP CUP NO:	Location of composite

NEW DISPLAY

NO OF STDS:	Total # of standards in run	
STD NO	CUP NO	The number under std no
1	:	refers to standard no..
2	:	Reply should indicate in
etc.		which cup std. no. 1 is
		located and so forth
STD DATA O.K., Y OR N:		Review the above table
		for accuracy. If changes
		are to be made type "N".
		Type "Y" to continue.
STD CONC:		Express in mg/ml
DIL FAC:		Dilution of tabs & comp (ml)
DECLARED IN MG OR GR?:		Type MG or GR (grains) to
		specify reporting units
AMT DECLARED:		Theo. amount/tablet, same
		unit as former.
WT OF COMP:		Composite weight (grams)
AVG TAB WR:		Express in grams
		(Note: If no composite
		type "1.0" in response
		to Wt of Comp and Avg
		Tab Wt.)

TAPE MOVEMENT

T:	Specify time interval
	in seconds at which the
	a/d is to be sampled.*

* A time interval of 0.3 secs. provides 2 - 3 points to define peak y maximum when data is collected from an auto-analyzer running 20 - 30 samples per hour.

6. After tapes stop and sample and display of data begins, depress sense switch 1 (SNS 1) to store initial baseline (BL1) and to begin search for peaks. (Note: Prior to depressing sense switch 1 data is not examined for peaks.)

7. The computer will now collect data and analyze it for peaks until SNS 1 is returned to its original position. This stops data collection, stores an average y value for the last two points to be used as BL2, and outputs results.

8. When calculations are complete the user is then asked "RECALC, Y OR N?". The user should now verify codes and determine if he desires to recalculate the results. If a "Y" is typed, %ONITA4 is loaded, the first absorbance value is outputted, and the user must type the proper code "s" (standard), "U" (unknown), "C" (composite) or "X" (deletion). Calculations are based on the user's coding. Deletion causes an absorbance to be omitted from the calculations.

9. After calculations are complete, the computer types the initial user input and asks the following questions. None of the following information is used by the program, but FOCAL-12 rules still apply to user input.

SPL NO:
DATE:
INITIALS:
PRODUCT:

(Note: In order to comply with FOCAL-12 rules governing the use of terminators within Ask statements. Mar. 12, 1972 may be written 12Mar72.)

10. After product is identified and the space bar is struck, the program stops and control returns to FOCALSNS. To restart the program type "L G, %ONITA1, Ø" and proceed as before.